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1. Lin Q Q, Zhou W Z, Wang Y F, et al. Abnormal peripheral neutrophil transcriptome in newly diagnosed type 2 diabetes patients[J]. *Journal of Diabetes Research*,2020: 1-10.
2. Sharma P, Singh S, Thakur V, et al. Novel and emerging therapeutic drug targets for management of type 2 Diabetes Mellitus[J]. *Obesity Medicine*,2021, 23:100329.
3. Dr Jayakrishnan B. Historical aspects and current understanding of the connections and implications of viruses and diabetes: A narrative review[J]. *Clinical Epidemiology and Global Health*,2022:1610 1110.
4. DiMeglio L A, Evans-Molina C, Oram R A. Type 1 diabetes[J]. *Lancet*,2018,391:2449-2462.
5. Bahmani M, Zargaran A, Rafieian-Kopaei Mahmoud. Ethnobotanical study of medicinal plants used in the management of diabetes mellitus in the Urmia, Northwest Iran[J]. *Asian Pac J Trop Med*, 2014, 7(Suppl 1): S348-S354.
6. Li Y G, Ji D F, Zhong S, et al. Hybrid of 1-deoxynojirimycin and polysaccharide from mulberry leaves treat diabetes mellitus by activating PDX-1/insulin-1 signaling pathway and regulating the expression of glucokinase, phosphoenolpyruvate carboxykinase and glucose-6-phosphatase in alloxan-induced diabetic mice[J]. *Journal of Ethnopharmacology*, 2011,134:961-970.
7. Dong J, Liang Q X, Niu Y. Effects of Nigella sativa seed polysaccharides on type 2 diabetic mice and gut microbiota[J]. *International Journal of Biological Macromolecules*, 2020,159 :725-738.
8. Castanheira F V S, Kubes P. Neutrophils and NETs in modulating acute and chronic inflammation[J]. *Blood*, 2019, 133(20): 2178-2185.
9. Soehnlein O, Steffens S, Hidalgo A, et al. Neutrophils as protagonists and targets in chronic inflammation[J]. *Nature Reviews Immunology*,2017, 17:248-261.
10. Garcia C, Fève B, Ferré P, et al. Diabetes and inflammation: fundamental aspects and clinical implications[J]. *Diabetes & Metabolism*,2010, 36:327-338.
11. Ilgun E, Akyurek O, Kalkan AO, et al. Neutrophil/Lymphocyte ratio and Platelet/Lymphocyte ratio in fibromyalgia[J]. *European Journal of General Medicine*,2016,13(2):100-104.
12. Wang Z, Gerstein M, Snyder M. RNA-Seq: a revolutionary tool for transcriptomics[J]. *Nature Reviews Genetics*,2009,10(1):57-63.

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